Bilirrubina

Standard

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Bilirrubina

Standard

To plot photocolorimetric calibration curves

PRINCIPLE

See the package insert for Wiener lab.'s **Bilirrubina**.

PROVIDED REAGENTS

S. Standard: 500 ug-vials of porcine ultrapure bilirubin according to AACC standards, to a final concentration of 100 mg/l.

NON-PROVIDED REAGENTS

- **Bilirrubina**, separately provided by Wiener lab.
- Distilled water.

INSTRUCTIONS FOR USE

- Remove the vial from its protective case. Take the metal precinct and open, slowly pulling out the rubber stopper to avoid any loss of material.
- Add exactly 5 ml distilled water.
- Cap, place in the case and mix by inversion (several times, for at least 30 minutes) until complete dissolution of the bilirubin.
- It is convenient to keep the Standard refrigerated during that process and afterwards.

WARNINGS

The Reagent is for "in vitro" diagnostic use.
Use the reagents according to the working procedures for clinical laboratories.
The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Bilirrubina Standard is stable at room temperature until the expiration date shown on the box. Once reconstituted, it is stable for 30 days in freezer or for 5 days in refrigerator (2-10°C). Protect it from light keeping it in its protective case.

INSTABILITY OR DETERIORATION OF REAGENTS

Any decoloration of the reagent indicates its deterioration. Moisture deteriorates the reagent.

REQUIRED MATERIAL (non-provided)

See the package insert for Wiener lab.'s **Bilirrubina**.

ASSAY CONDITIONS

See the package insert for Wiener lab's **Bilirrubina**.

PROCEDURE

Use Wiener lab.'s **Bilirrubina** reagents.
In 6 tubes labeled (1, B₁, 2, B₂, 3, B₃) place:

Tube	Standard	Reagent A	Water	Diazoreagent	Reagent B	Concentration (mg/l)
1	50 ul	2.6 ml	-	0.2 ml	-	25
B ₁	50 ul	-	2.6 ml	-	0.2 ml	-
2	100 ul	2.6 ml	-	0.2 ml	-	50
B ₂	100 ul	-	2.6 ml	-	0.2 ml	-
3	200 ul	2.5 ml	-	0.2 ml	-	100
B ₃	200 ul	-	2.5 ml	-	0.2 ml	-

Immediately mix each tube by inversion, after adding the Diazo reagent. After 5 minutes, read in spectrophotometer at 530 nm or in photocolorimeter with green filter at 520-550 nm, setting the instrument to zero with distilled water.

CALCULATIONS

Subtract to the readings of tubes 1, 2 and 3 their corresponding Blanks (B₁, B₂, B₃, respectively) obtaining the "corrected readings". With these values draw a graphic, in the usual way, having the corrected reading values and in the vertical axis and the bilirubin concentrations in mg/l on the horizontal axis. The reaction follows the Beer law, which allows to calculate a colorimetric factor for each tube.

$$\text{Factor} = \frac{\text{Bilirubin (mg/l)}}{\text{Corrected reading}}$$

Calculate the average factor.

PROCEDURE LIMITATIONS

Mistakes in the reconstitution may cause erroneous results.
See Procedure Limitations in the package insert of the kit in use.

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Wiener Laboratorios S.A.I.C.