Fer-color Transferrina

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- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

Fer-color Transferrina

Colorimetric method for the determination of serum
Total Iron Binding Capacity (TIBC)

SUMMARY

In the human body, iron circulates as Fe (III) bound to a specific carrier protein: transferrin or siderophilin. Its function is to bind iron from absorption sites (intestinal mucosa) or deposit (reticuloendothelial system) and to carry it to hematopoietic organs where it is used.

In a normal subject only a third part of transferrin is iron-saturated, the rest being free to bind and transport any eventual intake.

Transferrin physiological activity can be accurately determined by measuring the total iron binding capacity (TIBC). TIBC will be increased in posthemorrhagic and ferropenic anemias in general, in liver failure and physiologically in the last months of pregnancy.

On the other hand, TIBC is decreased in hemochromatosis, in certain types of anemias with disproteinemia (infectious, neoplastic, nephropathic, etc.) in chronic liver diseases and in great protein losses associated with nephrotic syndrome.

PRINCIPLE

Transferrin or specific iron carrier protein is assayed through its physiologic activity of binding Fe (III) (TIBC) at a pH higher than 7.2 in which transferrin is saturated in the presence of excess Fe (III). The remaining unbound Fe (III) is totally removed by coprecipitation with magnesium carbonate.

The iron bound to transferrin is released and colorimetrically determined according to **Fer-color** or **Fer-color AA** procedure. The quantity of transferrin is expressed as the micrograms of Fe (III) with which it is saturated.

PROVIDED REAGENTS

A. Reagent A: stabilized solution of Fe (III).

B. Reagent B: granulated magnesium carbonate.

NON-PROVIDED REAGENTS

- Wiener lab's **Fer-color** or **Fer-color AA** reagents.

- Distilled water.

INSTRUCTIONS FOR USE

Reagent A: ready to use.

Reagent B ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents are stable at room temperature until the expiration date shown on the box.

Reagent A: the container should be immediately capped after used.

INSTABILITY OR DETERIORATION OF REAGENTS

See Wiener lab's **Fer-color** or **Fer-color AA** manual instructions.

SAMPLE

See Wiener lab's **Fer-color** or **Fer-color AA** manual instructions.

REQUIRED MATERIAL (non-provided)

- See Wiener lab's **Fer-color** or **Fer-color AA** manual instructions.
- Kahn tubes.

ASSAY CONDITIONS

See Wiener lab's **Fer-color** or **Fer-color AA** manual instructions.

PROCEDURE

a) Transferrin saturation: in a Kahn tube place 500 ul serum and 500 ul Reagent A. Mix and let stand for 5 minutes at 37°C. Add the content of an Reagent B spoon. Cap and shake vigorously 5 minutes at room temperature. Shaking should be longitudinal. Centrifuge 10-15 minutes at 3,000-4,000 r.p.m. until a clear (or with serum own opalescence) supernatant is obtained.

b) Color reaction: procedure as indicated under Wiener lab's **Fer-color** or **Fer-color AA** manual instructions.

STABILITY OF FINAL REACTION

See Wiener lab's **Fer-color** or **Fer-color AA** manual instructions.

CALCULATIONS

Correct readings and perform calculations as indicated for serum iron determination, multiplying by two the final result due to serum dilution. Usually serum iron determination is performed along with transferrin determination. In that case three values are reported: Serum Iron, Transferrin (TIBC) and Transferrin Percent Saturation which is calculated as follows:

$$\% \text{ Saturation} = \frac{\text{Serum iron (ug/dl)}}{\text{Transferrin (TIBC) (ug/dl)}} \times 100$$

If serum iron determination is intended to be simultaneously performed so as to obtain the Percent Saturation Calculation, Transferrin Saturation should be started 30 minutes earlier.

REFERENCE VALUES

Literature report the following ranges for Transferrin and Transferrin Saturation in normal subjects:

Transferrin (TIBC): 250-400 ug/dl

Transferrin Saturation: 20-55 %

It is recommended that each laboratory establishes its own reference values.

PROCEDURE LIMITATIONS

See Wiener lab.'s **Fer-color** or **Fer-color AA** manual instructions.

Insufficient shaking will yield a false increase of transferrin values.

PERFORMANCE

a) Reproducibility: by replicates of the same sample in the same day the following data were obtained:

Level	S.D.	C.V.
280 ug/dl	± 23.2 ug/dl	8.29 %
560 ug/dl	± 28.9 ug/dl	5.16 %

b) Linearity: reaction is linear up to 500 ug/dl.

c) Detection limit: depends on the photometer used and the wavelength. According to the required sensitivity, for a minimum ΔA of 0.001, the minimum detectable concentration will be 1 ug/dl.

WIENER LAB PROVIDES

Auxiliary reagents to process 25 samples (TIBC):

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- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACCPress, 4th ed., 2001.

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Wiener Laboratorios S.A.I.C.