HDL Colesterol

Reactivo Precipitante

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HDL Cholesterol

Reactivo Precipitante

For HDL-cholesterol separation in serum or plasma

SUMMARY

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. These particles solubilize and transport cholesterol into the bloodstream. The relative proportion of protein and lipid determines the density of these lipoproteins.

The main role of high density lipoprotein (HDL) in lipid metabolism is the uptake and transport of cholesterol from peripheral tissues to the liver through a process known as reverse cholesterol transport (cardioprotective mechanism). Low HDL-cholesterol levels are associated with an increased risk of coronary heart disease. Therefore, the determination of serum HDL-cholesterol is a useful tool for identifying high-risk patients.

PRINCIPLE

High-density lipoproteins (HDL) are separated selectively precipitating the low and very-low-density lipoproteins (LDL and VLDL) by adding MW 50,000 dextran sulfate in the presence of Mg^{+} ions.

In the supernatant, separated by centrifugation, the HDL remain and the determination of cholesterol associated to them is performed using the enzymatic system cholesterol-oxidase/peroxidase with Trinder colorimetry (Phenol/4-Aminophe-nazone).

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INSTRUCTIONS FOR USE

Precipitating Reagent (A+B); preparation: in the provided bottle, measure 2.5 ml of Reagent A and 2.5 ml Reagent B. Mix by inversion, label and date.

Small quantities can be prepared according to what is needed, following the 1+1 proportion for both reagents.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at room temperature until the expiration date shown on the box.

Precipitating Reagent: is stable for 6 months at room temperature and 1 year in refrigerator (2-10°C) since preparation date.

INSTABILITY OR DETERIORATION OF REAGENTS

Any trace of bacterial contamination may be a sign of deterioration of the reagents.

SAMPLE

Serum or plasma

a) Collection: obtain the sample in the usual way.

b) Additives: if plasma is used, collect only with heparin.

c) Known interfering substances: anticoagulants other than heparin and bilirubinemia over 50 mg/l interfere.

See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: separate serum within one hour from extraction. The Lipid Research Clinics recommend refrigerating the sample until the test is performed. The storage or conservation of the samples at room temperature alters their lipoprotein composition, within 24 hours. Some authors mention a stability of 3 days at 4°C that is extended when frozen, but there is a large variability between different samples, therefore, it is recommended to keep the sample refrigerated and to process it within 24 hours.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer or photocolorimeter.
- Micropipettes and pipettes
- Kahn tubes.
- Tubes or spectrophotometric cuvettes.
- Water bath at 37°C.
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 505 nm in spectrophotometer or 490-530 nm in photocolorimeter with green filter.
- Reaction temperature: 37°C
- Reaction time: 45 minutes
- Sample volume: 500 μ l
- Precipitating Reagent volume: 50 μ l
- Supernatant volume: 100 μ l
- Working Reagent volume of **Colestat enzimático** or **Colestat enzimático AA/Línea Líquida**: 2 ml
- Final reaction volume: 2.1 ml

PROCEDURE

In a Kahn tube measure 0.5 ml (500 ul) sample and add 50 ul Precipitating Reagent. Homogenize by shaking (not by inversion) for 20 seconds and place for 30-40 minutes in refrigerator (2-10°C) or 15 minutes in water bath at that temperature. Do not freeze. Centrifuge 15 minutes at 3000 rpm. Use the clean supernatant as sample.

In three tubes labeled B (Blank), S (Standard) and U (Unknown) place:

	B	S	U
Supernatant	-	-	100 ul
Standard	-	20 ul	-
Working Reagent	2 ml	2 ml	2 ml

Mix and incubate for 5 minutes at 37°C, if the Working Reagent of **Colestat enzimático AA/Línea Líquida** is used or 15 minutes at 37°C if **Colestat enzimático** Working Reagent is employed. Remove from bath and let cool. Read at 505 nm in spectrophotometer or at 490-530 nm in photocolimeter with green filter, setting the instrument to zero with the Blank.

STABILITY OF FINAL REACTION

The reaction color is stable for 2 hours, thus the absorbance must be read within that period.

CALCULATIONS

$$\text{HDL Cholesterol (g/l)} = U \times f \quad f = \frac{0.457}{S}$$

$$0.457 = 2 \text{ (g/l)} \times \frac{FV_E}{V_S} \times \frac{RV_E}{RV_{Std}} \times \frac{V_{Std}}{V_E}$$

FV_E = extract final volume = 0.55 ml

V_S = processed sample volume = 0.5 ml

RV_E = reaction volume with extract = 2.1 ml

RV_{Std} = reaction volume with Standard = 2.02 ml

V_{Std} = Standard volume in the reaction = 0.020 ml

V_E = extract volume in the reaction = 0.1 ml

If Reagent volumes other than 2 ml are used, the 0.457 factor varies and must be calculated again, replacing VR_E and VR_S in the formula.

REFERENCE VALUES

The National Cholesterol Education Program (NCEP) Expert Panel provided the following reference values:

0.40 - 0.60 g/l

It is recommended that each laboratory establishes its own reference values. However, values over 0.40 g/l are recommended and the ones over 0.60 g/l have been regarded as protective. On the contrary, HDL-cholesterol values below 0.40 g/l are considered as an indicator of coronary heart disease risk.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

The accuracy and precision of the determination mainly depend on following the precipitation conditions, thus, the established times and temperatures must be followed, although they do not need a strict control.

When the HDL cholesterol cannot be completely separated in a common centrifuge, due to high levels of triglycerides, the determination can be performed as follows:

Follow the instructions under PROCEDURE up to the incubation at 4-10°C, then place the reaction mixture in capillary tubes and centrifuge in a microhematocrit centrifuge at 10,000 rpm for 5 minutes. Cut the tube discarding the precipitate and use the clean supernatant for the test.

PERFORMANCE

a) Reproducibility: processing replicates of one sample on the same day, the following results were obtained:

Level	S.D.	C.V.
0.29 g/l	± 0.011 g/l	3.8 %
0.63 g/l	± 0.023 g/l	3.7 %

b) Linearity: the reaction is linear up to 5 g/l.

c) Detection limit: depends on the spectrophotometer used. For a 0.001 O.D. reading, the minimum visible change of concentration will be of approximately 0.0063 g/l.

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Wiener Laboratorios S.A.I.C.