Bilirrubina

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Bilirrubina

For direct and total bilirubin determination

SUMMARY

Bilirubin, compound produced by degradation of hemoglobin, is captured by the liver for its conjugation and bile excretion. Liver cell disorders or bile obstructions may lead to hyperbilirubinemia. *Erythroblastosis fetalis* or hemolytic anemia of the newborn is a pathology caused by mother-fetus blood incompatibility, producing an excessive destruction of red blood cells. This results in a severe increase in serum Bilirubin with risk of pigment diffusion to the central nervous system. Consequently, bilirubin determination in these newborns is highly important.

PRINCIPLE

Bilirubin reacts specifically with diazotized sulfanilic acid producing a purple-red pigment (azobilirubin) photocolormetrically measured at 530 nm. While conjugated (direct) bilirubin reacts directly with the diazoreagent, unconjugated (indirect) bilirubin requires the presence of an aqueous developer (Reagent A) to enable reaction. Therefore, caffeine benzoate must be added to the reaction medium so that the total bilirubin (conjugated and unconjugated) in the sample generates reaction.

PROVIDED REAGENTS

- A. Reagent A:** 0.13 mol/l caffeine benzoate aqueous solution, buffered and stabilized.
B. Reagent B: 29 mmol/l sulfanilic acid solution and 0.17 mol/l hydrochloric acid.
C. Reagent C: 0.07 mol/l sodium nitrite solution.

NON-PROVIDED REAGENTS

- Distilled water.
- Wiener lab's **Bilirrubina Standard** for periodic instrument calibration.

INSTRUCTIONS FOR USE

Reagent A: ready to use.

Reagent B: ready to use.

Diazoreagent: according to the volume to prepare, mix 1 part of Reagent C with 21 parts of Reagent B. Label and date.

WARNINGS

The reagents are for «in vitro» diagnostic use. Use the reagents maintaining the usual working precautions for the clinical laboratory. All reagents and samples must be discarded according to local regulations in force.

Reagent B: H315+H320: Causes skin and eye irritation. H314: Causes severe skin burns and eye damage. P262: Do not get in eyes, on skin, or on clothing. P305 + P351 +

P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P280: Wear protective gloves/protective clothing/eye protection/face protection.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at room temperature (< 25°C) until the expiration date shown on the box.

Diazoreagent: at 2-10°C and in caramel-colored glass bottle it is stable for 3 months from the date of preparation.

INSTABILITY OR DETERIORATION OF REAGENTS

The presence of sediment and/or discoloration of the reagents may indicate reagent deterioration.

SAMPLE

Serum, plasma or amniotic fluid

a) Collection: obtain serum or plasma as usual. Protect from natural or artificial light covering the tube with black paper. You can also make the determination in amniotic fluid.

b) Additives: use heparin for collection in case the sample to be used is plasma

c) Known Interfering Substances: samples with mild or intense hemolysis inhibit the direct reaction, obtaining falsely increased total bilirubin values. Refer to Young, D.S. in references for the effect of drugs on the present method.

d) Stability and storage instructions: sample should be preferably fresh. If the assay is not performed immediately, serum can be stored up to 48 hours in refrigerator (2-10°C) and whole blood no more than 24 hours in refrigerator (2-10°C) or 12 hours at room temperature. Keep the amniotic fluid frozen until the assay is performed.

The action of light is capable of destroying up to a 50% of the bilirubin present in the sample in one hour. Consequently, it should be carefully protected.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer and photocolormeter
- Micropipettes and pipettes for measuring stated volumes
- Caramel-colored glass bottle
- Tubes
- Watch or timer

ASSAY CONDITIONS

- Wavelength: 530 nm in spectrophotometer or 520-550 nm in photocolormeter with green filter
- Reaction temperature: room temperature

- Reaction time: 5 minutes
- Sample volume: 200 ul
- Final reaction volume: 2.9 ml

PROCEDURE

I-TECHNIQUE FOR SERUM

In three tubes labeled B (Blank), D (Direct) and T (Total) place:

	B	D	T
Sample (serum)	200 ul	200 ul	200 ul
Distilled water	2.5 ml	2.5 ml	-
Reagent A	-	-	2.5 ml
Reagent B	200 ul	-	-
Diazo reagent	-	200 ul	200 ul

Mix immediately by inverting each tube. After 5 minutes, read on spectrophotometer at 530 nm or on photo-colorimeter with green filter (520-550 nm), taking the instrument to zero with distilled water. Readings can be performed at intervals between 4 and 15 minutes, except for direct bilirubin, which must be read at exactly 5 minutes. If read beforehand, there will be undervaluation of the results by incomplete reaction. If read afterwards, there will be overvaluation because free bilirubin starts reaction. Icteric sera: the technique described should be used with smaller amounts of sample, according to the severity of the jaundice. Thus, for moderate jaundice 50 ul serum will be used, while for severe jaundice only 20 ul will be required. Multiply the results by 3.79 and 9.38 respectively.

II- TECHNIQUE FOR AMNIOTIC FLUID

Total bilirubin is determined. In two tubes labeled B (Blank) and T (Total) place:

	B	T
Sample (amniotic fluid)	1 ml	1 ml
Distilled water	1.5 ml	-
Reagent A	-	1.5 ml
Reagent B	0.2 ml	-
Diazo reagent	-	0.2 ml

Proceed in the same way as in Technique I. Divide the result by 5.37.

CALCULATIONS

Total Bilirubin (mg/l) = (T - B) x f

Direct Bilirubin (mg/l) = (D - B) x f

Free Bilirubin (indirect) = Total BRB - Direct BRB

The colorimetric factor (f) must be calculated using Wiener lab's **Bilirubina Standard**.

QUALITY CONTROL METHOD

If the test sample is serum or plasma, process 2 levels of a quality control material (**Standatrol S-E 2 niveles**) with known concentrations of bilirubin, with each test.

REFERENCE VALUES

Bilirubin in serum or plasma

- Adults:

Direct: up to 2 mg/l

Total: up to 10 mg/l

- Newborns:

	Full-term infant	Premature
Cord blood	25 mg/l	
Up to 24 hrs	60 mg/l	80 mg/l
Up to 48 hrs	75 mg/l	120 mg/l
3rd to 5th day	120 mg/l	240 mg/l

The values begin to decrease to reach the average adult level after one month of birth. In premature infants, bilirubin levels take longer to reach normality, depending on the degree of hepatic immaturity.

Bilirubin in amniotic fluid

Values lower than 1 mg/l are considered normal, while between 1 and 2.7 mg/l suggest the possibility of an affected fetus.

Levels greater than 2.7 mg/l are found only in presence of fetal erythroblastosis. If the amount of bilirubin is greater than 4.7 mg/l the fetus is affected and probably has some kind of circulatory insufficiency. If the amount of bilirubin exceeds 9.5 mg/l fetal death is imminent.

It is recommended that each laboratory establish its own reference values.

UNITS CONVERSION

Bilirubin (mg/l) = Bilirubin (mg/dl) x 10

Bilirubin (umol/l) = Bilirubin (mg/dl) x 17.1

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE. The action of light, both on sera as on the standard solutions, is capable of destroying up to 50% of the existent bilirubin in one hour.

PERFORMANCE

a) Reproducibility: processing replicates of the same samples in the same day, the following values were obtained:

Direct Bilirubin

Level	S.D.	C.V.
2.1 mg/l	± 0.18 mg/l	8.5 %
26.0 mg/l	± 0.62 mg/l	2.4 %

Total Bilirubin

Level	S.D.	C.V.
9.8 mg/l	± 0.47 mg/l	4.7 %
127.7 mg/l	± 1.45 mg/l	1.1 %

b) Linearity: reaction is linear up to 150 mg/l.

c) Recovery: by adding known amounts of bilirubin to various

sera, a recovery between 100% and 106% was obtained.
d) Detection limit: it depends on the spectrophotometer used and the wavelength. According to the sensitivity required for a minimum ΔA of 0.001, the smallest detectable concentration change will be 0.17 mg/l.

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Wiener Laboratorios S.A.I.C.